



BIOSYNERGY, INC.

1940 E. Devon Avenue, Elk Grove Village, Illinois 60007
(847) 956-0471 Toll Free (800) 255-5274 Fax (847) 956-6050
e-mail: biosynergy@att.net website: www.biosynergyinc.com

HEMOTEMP® II ACTIVATOR

Biosynergy, Inc. introduces its HemoTemp® II Activator, a block model heater designed to hold and heat as many as 12 HemoTemp® II indicators. The Activator is designed to be a companion to the HemoTemp® II indicators. Since the HemoTemp® II indicator's irreversible portion needs to be heated immediately before placing on a cold blood bag to be fully time/temperature activated, the Company designed the Activator as a precise and accurate heat source for this portion of activation.

The HemoTemp® II Activator eliminates the uncertainty about whether the irreversible portion of the indicator has been properly heat activated. The system offers the flexibility of storing up to 12 indicators that can remain in the unit for as long as 4 weeks. The unit is portable and simple to operate, just plug in, and turn on.

The HemoTemp® II indicators are simply placed with the irreversible ("window") portion down into the Activator thus heating the indicator to the precise temperature necessary for proper heat activation. The Activator offers temperature uniformity of +/- 0.5°C with a temperature stability of +/- 2°C. The unit is **UL Approved** and comes with a **2 year Manufacturers' Warranty**.

The HemoTemp® II Activator is ideal for new HemoTemp® II indicator users or those who have been using the HemoTemp® II indicator for years. If you would like additional information, please feel free to call our toll-free number, **800-255-5274**, or email biosynergy@att.net and we would be happy to assist you.