

**CONTROL ENZYME LOSS IN BLOOD SAMPLES**


Use **LabTemp®20** To Control Temperature Base

2	BSI
5	20
8	17
11	14


LabTemp® 20 Indicator

2.0°C  
(2 is GREEN)

LabTemp® 20



**EFFECT OF TEMPERATURE**



Typical progress curves of an enzyme reaction at various temperatures, showing dependence of apparent optimum temperature on time.

**COLOR INTERPRETATION**

- **TAN/RED:**  
Temperature is 1°C lower than the digit that appears
- **GREEN:**  
Temperature is indicated by the digit that appears
- **BLUE:**  
Temperature is 1°C higher than the digit that appears

